




DATA SHEET 3942

VORTRAN APM-Plus Connection Kit (VCK-3942) For Use with VORTRAN APM-Plus



Description:

Connection Kit 3942 is intended for use with the VORTRAN APM-Plus. The Connection Kit connects to a breathing circuit using an included tee adapter (male and female 22mm conical connector as per ISO 5356-1). The VORTRAN APM-Plus is used to measure airway pressure for any patient utilizing positive pressure devices where monitoring is desired.

Dimensions and Weight:

1.09 ounces, 31 grams
8" x 10" packaged. Included tubing is 4 Feet.

Application:

Single-patient-use only. Disposable. Do not reuse. Disassembly, cleaning, sterilization, maintenance, and factory service are not required. Max use 30 days.

Compatible with:

APM-Plus 3960

Components / Materials:

Luer Adapter (PN 3911) / Nylon
Bacteria Filter (PN 3910) / Housing: Acrylic and Cyrolite G20-Hiflo, Filter Media: PTRE
Oxygen Tubing, 4' (PN 3909) / Polyurethane
Tee Adaptor (PN 3933) / Polypropylene
All materials are DEHP and natural latex free.

Storage:

Temperature range of -25° to 70° Celsius
Humidity range of 15 to 90%

Disinfection and Sterilization:

Once used, do not reuse, clean, or sterilize it. Single-patient-use and disposable. Packaged non-sterile.